

# Data Quality Monitoring for High Energy Physics (DQM4HEP) Module interfaces

**R. Été, A. Pingault, L. Mirabito**

Université Claude Bernard Lyon 1 - Institut de Physique Nucléaire de Lyon / Ghent University

10 février 2016




Université Claude Bernard



Lyon 1



## Global workflow




# Module applications - analysis module

## Purpose

- Receive events from a collector server and process them
  - Produce monitor elements (histograms, scalars, generic TObject)
  - Follow the run control signals (SOR, EOR)
- 
- **Init** : Initialize the application : load dlls, declare services, etc ... Wait for a SOR
  - **Start of run** : start cycles loop, open archive
  - **Start of cycle** : start a cycle of '*process event*'
  - **Process event** : Process incoming event, fill monitor elements, etc ...
  - **End of cycle** : send subscribed monitor elements, update archive (opt).
  - **End of run** : Wait for SOR, close archive (opt).
  - **End** : Clean and exit module.


To implement online DQM analysis, user must implement the `DQMAnalysisModule` interface. A shared library must be build and loaded in the application using the plugin system (see next slides).

Use `dqm4hep_start_analysis_module` to start an analysis module.



Analysis module  
application flow

## Module API



# Module applications - standalone module


## Purpose

- No event reception
  - No run signals
  - Produce monitor elements (histograms, scalars, generic TObject)
- 
- **Init** : Load dlls, init the module.
  - **Start of cycle** : start a timer cycle of n seconds
  - **Process** : call back function.
  - **End of cycle** : collect monitor elements and send
  - **End** : The application has received a signal to exit and the process ends.


To implement online standalone analysis, user must implement the `DQMStandaloneModule` interface. A shared library must be build and loaded in the application using the plugin system (see next slides).

Designed for *slow control - like* data processing.

Use `dqm4hep_start_standalone_module` to start a standalone module.



Standalone module application flow



Data processing performed in **modules** (standalone or analysis).

Modules **book** monitor elements, **fill** them and **publish** them to a single collector.

A monitor element is a wrapper around a ROOT TObject with some additional attributes :

→ Type, name, path (i.e "/Efficiency/Layer2/"), collector name, quality flag, reset policy, title, description, run number, quality test results.

The `DQMModuleApi` class provide a static interface to perform operations within the application :

- Monitor elements management (book, delete, reset, from xml)
- Directory structure management (mkdir, cd, ls, rmdir, pwd)
- Quality test management (register, add, remove, run, from xml)

Quality tests can be run on a particular monitor element to test the quality of the processed data (chi2, Kolmogorov, user defined).

Note that **QTest results are sent to the collector together with the monitor element !**

## Example modules

An LCIO example module can be found on the github page of DQMCore package

<https://github.com/DQM4HEP/DQMCore> :

CaloHitModule (lcio, analysis) :

- `source/examples/module/lcio/CaloHitModule.h`
- `source/examples/module/lcio/CaloHitModule.cc`
- `conf/lcioCaloHit.xml`
- `dqm4hep_start_analysis_module` **executable**